# Venus® Bulk Flow ONE

Venus® Bulk Flow ONE

**DE** Gebrauchsanweisung  
**GB** Instructions for use  
**FR** Mode d'emploi  
**ES** Instrucciones de uso  
**IT** Istruzioni per l'uso  
**PT** Instruções de uso

**NL** Gebruiksaanwijzing  
**SE** Bruksanvisning  
**DK** Brugervejledning  
**NO** Bruksanvisning  
**FI** Käyttöohjeet

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ONE 20s 4mm

7. Alternativ, kann Venus Bulk Flow ONE can applied in a thin layer as a liner on exposed dentine or used to repair direct and indirect restorations in combination with a suitable adhesive (e.g. BOND® Adhesive Systems, see manufacturer's instructions for use to apply).

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12. At the end of treatment, remove the hygienic protective cover from the syringe and discard it.

13. Do not reuse the syringes if there is evidence of contamination with body fluids.

14. Discard the used PLT.

**Tips on product use**

Venus Bulk Flow ONE should be processed at an ambient temperature of 23°C (73°F).

1. The polymerisation times indicated above are minimum times. They should not be reduced under any circumstances to ensure that an optimum polymerisation rate and the associated physical properties of the material thereof are achieved.

2. When building up large restorations, or those which require a lengthy amount of time to complete, the overhead operating light intensity must be turned off or moved away from the working area in order to prevent premature curing of the composite.

3. Deep cavities where there is only a minimal layer of dentine over the pulp or where a small area of the pulp is already exposed must be capped locally with a calcium hydroxide product either indirectly or directly. The calcium hydroxide product should be capped locally with a glass ionomer cement. Venus Bulk Flow ONE is placed in the cavity. Prepared proximal surfaces and occlusal margins are smoothed with a dental bur. The cavity is then rinsed with water and dried with a dry cloth. The cavity is then filled with the composite resin. Avoid contact with products containing eugenol, as they may limit the optimum polymerisation of the composite.

4. When using Venus Bulk Flow ONE PLTs, use a gentle, even motion when exerting pressure on the applicator gun. Use of excessive force or sudden movement could result in unanticipated extrusion of the material or ejection of the PLT from the applicator gun, creating a potential hazard.

5. It is also possible to apply the Venus Bulk Flow ONE in fine couches liner on the dentine surface or used to repair direct and indirect restorations in combination with a suitable adhesive (e.g. BOND® Adhesive Systems, see manufacturer's instructions for use to apply).

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9. Bitte kann Venus Bulk Flow ONE can applied in a thin layer as a liner on exposed dentine or used to repair direct and indirect restorations in combination with a suitable adhesive (e.g. BOND® Adhesive Systems, see manufacturer's instructions for use to apply).

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12. At the end of treatment, remove the hygienic protective cover from the syringe and discard it.

13. Do not reuse the syringes if there is evidence of contamination with body fluids.

14. Discard the used PLT.

**General information**

Please refer to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

**Warnings and Safety instructions**

May cause skin and eye irritation. Avoid contact with skin. Skin contact may cause sensitization. Wear suitable protective gloves in case of contact with skin. Immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist. May damage fertility.

For SSCP please see the EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact sscp-service@kulzer-dental.com. To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product.

Stand: 2023-09

**■ Manufacturer:**  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)  
Made in Germany

Distributed in USA/Canada exclusively by:  
Kulzer, LLC  
4315 South Lafayette Blvd.  
South Bend, IN 46614-2517  
1-800-431-1785

**Caution:**  
Federal law restricts this device  
to sale by or on the order of a dentist.

**Dentistry – Polymer-based restorative materials**  
ISO 4049 (type 1, class 2, group 1)

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