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STERILE R

C-TECH IMPLANT s.r.l.  
Via Cesare Battisti n. 2,  
40123 Bologna (BO) - Italy

**INSTRUCTIONS FOR USE  
POSITIONING OF THE IMPLANTS**

**DESCRIPTION** - The C-TECH IMPLANT "GL", "BL", "ND" and "EL" dental implant systems are made for use in a total or partial, mandibular or maxillary edentulism. It is used as a support for a fixed or removable prosthesis.

**SELECTION OF THE PATIENT** - The assessment of the patient before a surgical implant procedure is extremely important. This includes the assessment of the patient's general state of health, state of oral hygiene, good dental and periodontal condition, and anatomical acceptability. A case history and clinical history of the patient is essential.

**PRE-OPERATIVE TREATMENT** - The choice of the size of the implant is important for long-term success. It is better to use the maximum diameter and length possible, for greater stability of the superimposed prosthesis. The measurements can be made directly on the panoramic x-ray, superimposing the transparencies on it, enlarged by 25% compared to the real image. The edges of the crest must be adequately palpated in order to estimate an insertion angle that must reach a good parallelism with the other implants, and have an emergence which allows for proper prosthetic rehabilitation.

**SURGICAL PROCEDURE** - As with any other surgical action, it is important for the implant procedure to be performed in an environment which is as sterile as possible! The surgical kit is not sterile when supplied and must be autoclaved before use. Follow the sterilization cycle instructions given in paragraph *Sterilization phases* in the instruction sheet F017

**PREPARATION OF THE SITE** - Make a full-thickness mesiodistal incision along the surface of the crest. Using a periosteal elevator, leave the surface exposed to the alveolar bone only as necessary in order to have adequate vision and adequate working space for the surgical procedure. Pointed crests or other bone irregularities must be removed using a rose bur or a bone forceps. A minimum distance of 4-6 mm (from cut to cut) must be maintained between the implants and/or the adjacent natural teeth.

**GENERAL INSTRUCTIONS FOR PREPARING THE IMPLANT SITE**

It is compulsory for the entire osteotomy cut procedure to be carried out at a speed of 10 to 1000 RPM. Abundant external irrigation is required for cooling.

**STEP BY STEP PREPARATION OF THE IMPLANT SITE** - The implant site must be prepared in a visible, clear field, so that the operator can prepare it with precision in every phase, thus allowing for proper positioning of the implant.

For the preparation of the implant site, it is necessary to use a series of burs with external cooling used sequentially. For every depth and widening drill (except Ø 5.1) are available removable stops in steel which can facilitate the implant site preparation. Phase 4 (preparation of the head) can be also done at the beginning of the surgery depending on the chosen surgical protocol.

- FIRST DRILL** (guide drill or lance drill) allows for the direct perforation of the bone cortex at the centre of the crest.) Operating speed 100-750 RPM
- SITE PREPARATION IN LENGTH** The depth drill, size 2.2, equipped with mar-

**DENTAL IMPLANT SYSTEMS**

ked stops indicating the depth, prepares the implant site of the selected length. Operating speed 100-750 RPM.

- SITE PREPARATION IN WIDTH** The shaping drills for the diameter of the selected implant are calibrated with the core of the screw. They do not cut at the tip, and are to be used in sequential manner with sizes 3.5, 4.0 to 5.0 in order to obtain the maximum precision with the minimum heating possible. For ND implant insertion use drill Ø 2.6 after drill Ø 2.2. Operating speed 100-750 RPM.
- PREPARATION OF THE HEAD** For Ø 3.5 and 4.3 the preparation is finalized with counterbore drill which will adapt the inclination related to the chosen implant. If needed, use counterbore Ø 3.0 for insertion of ND implants. Operating speed 100-750 RPM.

**TREATMENT AND STERILIZATION** - The implants are supplied in sterile packages. Given the nature of surfaces treatment (sandblasting and acidification) of C-TECH IMPLANT "GL", "BL", "ND" and "EL" special care is required in handling it. Only neutral gloves or titanium instruments must be used to handle C-Tech Implants.

- Sterilization: the surgical instruments are not sterile, and must be sterilized before use.

**POSITIONING OF THE IMPLANTS**

- Irrigate the entire implant site with sterile water or saline solution.
- Bring the implant to the already prepared site. The implant must be brought to the sur-gery site without being removed from its cap. Great care is required in handling the implant.
- Using your fingers to put pressure on the cap, partially position the implant in the site. Secure the implant with the contrangle or the manual one, then potentially with the ratchet wrench, rotating the plant clockwise until it is fixed. The implant is properly position-ed when the platform is 1 mm above the crest of the alveolar bone. A cap screw is included above the cap of the implant, and this must be tightly screwed onto the implant before suturing the flap.
- An x-ray check of the position of the implant is recommended before closing the flap.
- The mucoperiosteal flap must be delicately repositioned for optimal adaptation of the tissue, and then sutured.

**POST-OPERATIVE COURSE** - The patient must be taught to follow the postoperative course, which includes cold compresses for the first 24 hours. An antibiotic must be prescribed an hour before the operation, and for four days after the operation. The stitches must be removed after one week. It is recommended that every removable prosthesis which remains near the implant site be adequately raised and aligned using a soft material per regeneration. Insufficient available bone, poor bone quality, poor oral hygiene on the part of the patient, tobacco smoking and systemic and general diseases (diabetes, rheumatoid arthritis, etc.) can contribute to poor bone integration and a consequent failure of the implant.

**SUPPLY** - C-TECH IMPLANT "GL", "BL", "ND" and "EL" dental implant systems implants are available in sterile packages with surface treatment. Each package contains an internal label. That label must be affixed to the patient's file for future reference. The label indicates the specifications of the positioned implant, including the date of reference. The following sizes are available:

**DIAMETER LENGTH**

- Ø 3.0 mm L. 9/11/13 mm
- Ø 3.5 mm L. 7/ 8/ 9/11/13/ 15 mm
- Ø 4.3 mm L. 6/ 7/ 8/ 9/ 11/ 13/ 15 mm
- Ø 5.1 mm L. 6/ 7/ 8/ 9/ 11/ 13/ 15 mm

**CONTRAINDICATIONS:**

Complications: lack of bone integration is shown by infections, movement or loss of bone. Each failed implant must be removed as soon as possible, and all of the surrounding tissue must be removed.

Caution: The surgical techniques required in order to position the dental implants are highly specialized and complex procedures. Therefore, a specialization course is highly recommended. The professionals must participate in and complete courses of study in order to be prepared for implant techniques; Improper techniques can cause the failure of implants, loss of bone, and postoperative complications.

Precautions: The adequacy of the bone must be determined by x-rays, palpations and by means of visual inspections of the site selected for the implant. Determine the location of all of the anatomic structures that must be avoided before beginning any type of implant procedure.

Warnings: Under law, only surgeons and certified dentists are allowed to use implants in surgery. Any use by other persons is strictly prohibited.

**CONTRAINDICATIONS**

**ABSOLUTE CONTRAINDICATIONS** - Known allergies to titanium, recent myocardial infarction, valvular prosthesis, severe renal disorder, severe liver disorders, osteomalacia, treatment-resistant diabetes, recent high-dosage radiation therapy, severe chronic alcoholism, drug abuse, illnesses and tumours in the terminal phase.

**GENERAL CONTRAINDICATIONS** - Non-compliant patient. Intraoral contraindications, such as: Anatomical impairments, insufficient bone tissue or poor-quality bone tissue.

**RELATIVE CONTRAINDICATIONS** - Chemotherapy, moderate renal disorders, moderate liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding or motivation, AIDS, HIV-seropositivity, prolonged cortico-steroid use, calcium-phosphorus metabolism disorders, erythropoietic disorders.

**DENTAL IMPLANT KIT**

Caution: this product is for single use only the device must be used exclusively by the physician. The "use by" expiration date is valid only if the product is stored properly in unopened, undamaged packaging. Disposable. Do not re-use and discard after use in accordance with the applicable regulations. The re-use involves risk of cross-infections and loss of device performances.

**ATTENTION PLEASE: C-TECH, declines every responsibility for any improper use of device unless instructed on enclosed directions leaflets or for use by nonprofessionals or not trained professionals, which may hurt or put under serious risks patient's health.**

**C-TECH IMPLANT s.r.l.**  
Via Cesare Battisti n. 2,  
40123 Bologna (BO) - Italy